

Your Partner in Laboratory Centrifuges



Micro Centaur Plus

USER MANUAL

PLEASE READ THIS BEFORE USE

MSE SAS

Health and Safety at Work

MSE SAS as designers, manufacturers, suppliers and importers of articles for use at work has to ensure that, as far as is reasonably practicable, articles which we design, produce, supply or import are safe and without risk to health.

We are required to provide information on the safety and handling precautions to be observed when installing, operating, maintaining and servicing our products. Such advice is contained in this manual.

We are also obliged to update this information should circumstances change and to operate a system to this end.

We should also like to point out, however that **you as users have an important responsibility in the provision and maintenance of safe working practices and conditions.**

Accordingly, we draw the following matters to your attention:

- 1. This apparatus should only be used as intended, within its design parameters by suitably qualified and trained personnel who have read and understood the relevant sections of this manual.**
- 2. This manual should be readily available to such personnel at all times.**
- 3. In addition to that which is written in the manual, normal common-sense safety precautions must be taken at all times to avoid the possibility of accidents. Particular care is required when working with apparatus at high temperature or pressure.**
- 4. Installation, maintenance, repairs and servicing should only be carried out by an MSE SAS approved engineer, and connection to electrical supplies should only be carried out by suitably trained personnel.**

TECHNICAL SUPPORT, WARRANTY SERVICE AND MAINTENANCE

If you are in any doubt whatsoever regarding the correct use of this apparatus, or if you require any technical data or assistance, please contact the MSE SAS FR Technical Support Department at:

MSE SAS

**2 rue des Artisans
49340 NUAILLE
FRANCE**

**Telephone: +332 41 64 73 41
E-mail: info@mse-sas.fr
Web Page: www.mse-sas.com**

OVERSEAS CUSTOMERS: Should contact their local MSE SAS Distributor, details can be found on our website.

ELECTRICITY SUPPLIES: Voltage and frequency

MSE SAS electrical apparatus is offered and labelled for one, or for a choice of two or more voltage ranges and, where necessary, different frequencies of mains supply. MSE SAS does not accept any responsibility for the operation of any such apparatus should it be connected to electricity supplies which are normally outside, or vary outside, the stated voltage and frequency values for which it is designed, nor for any consequential loss, damage or injury, howsoever caused.

Thank you for buying a **MSE MICROCENTAUR PLUS Centrifuge**. Please read this operating manual before using your centrifuge; it will provide you and your colleagues with useful information on all aspects of the centrifuge and its accessories.

MSE SAS may change the manual at any time and without notice because of improvements, typographical errors or improvements to facilities

All rights reserved. No part of this User Manual may be modified, distributed, published or reproduced without the prior permission of MSE.

You can find the current version of the user manual on our website under: www.mse-sas.com

As our customer, we should like to ensure that you are totally satisfied at all times.

Do not hesitate to contact us.

Your views are very important to us.

Warranty

The warranty period is **24 months** (unless otherwise specified in the purchase documents).
The service life of the centrifuge specified by the manufacturer is **10 years**.

After termination of the warranty period, it is necessary to carry out yearly technical inspections of the centrifuge.

The Manufacturer reserves the right to make technical changes in manufactured products.

The maximum period of storage of for centrifuges that are not used is **1 year**. After this period, a technical inspection of the centrifuge should be carried out by service personnel authorized by the manufacturer.

Visit our website : www.mse-sas.com for registration.

Warning signs:





	WARNING! Warning of potential injury or health risk.
	DANGER! Risk of electric shock with potential for severe injury or death as a consequence.
	DANGER! Biohazard with potential for risk to health or death as a consequence.
	DANGER! Risk of explosion with potential for severe injury or death as a consequence.

TABLE OF CONTENTS

1. APPLICATION	6
2. TECHNICAL DATA	7
3. INSTALLATION.....	8
3.1 UNPACKING THE CENTRIFUGE.....	8
3.2 LOCATION	8
3.3 CONNECTION TO THE POWER SUPPLY.....	8
3.4 FUSES	8
4. SAFETY NOTES.....	9
4.1 GENERAL INFORMATION ABOUT YOUR CENTRIFUGE	9
4.2 OPERATING PERSONNEL.....	9
4.3 SAFETY DEVICES.....	9
4.3.1 Lid lock.....	9
4.3.2 Checking the rotor has stopped.....	9
4.4 OPENING THE LID MANUALLY	9
4.5 GOOD OPERATOR PRACTICE	10
4.6 HAZARDS AND PRECAUTIONS	10
4.7 STORAGE	11
4.8 MAXIMUM LOAD.....	11
5. OPERATING THE CENTRIFUGE	12
5.1 GETTING STARTED	12
5.2 CONTROLS	12
5.2.1 The control panel.....	12
5.3 DESCRIPTION OF SOUNDS.....	14
5.4 LCD DISPLAY	14
6. MAINTENANCE.....	17
6.1 GENERAL CLEANING OF THE CENTRIFUGE	17
6.2 CLEANING OF THE ACCESSORIES	17
6.3 STERILIZATION AND DISINFECTION OF THE BOWL AND ACCESSORIES.....	17
6.4 TUBE CRACKING	17
7. SERVICING YOUR CENTRIFUGE.....	18
7.1 FAULT FINDING	18
7.2 OPENING THE LID MANUALLY	18
7.3 SAFETY INSPECTIONS	18
7.4 NOTES FOR ENGINEERS	18
8. HOW TO OBTAIN SERVICE ON YOUR MICRO CENTAUR PLUS.....	19
8.1 WARRANTY	19
8.2 MANUFACTURER'S INFO	19
9. ACCESORIES LIST	1
MICRO CENTAUR PLUS SERVICE REQUEST FORM	3
EU DECLARATION OF CONFORMITY	4
NOMOGRAPH.....	5

1. APPLICATION

The **MSE MicroCentaur plus** is a laboratory micro centrifuge intended for *in vitro* diagnostic (IVD). Its construction ensures easy and safe operation, for a wide range of applications in laboratories engaged in routine medical analysis and biochemical research works.

It is intended for separation of mixtures, suspensions, and systemic fluids into constituents of different densities under influence of centrifugal force. This centrifuge is an air-cooled, brush-less instrument and therefore during centrifugation the samples should be spun in sealed containers and sealed rotors.

This centrifuge is not intended to be used to centrifuge caustic, inflammable and explosive preparations.

2. TECHNICAL DATA

Type	MSB010.CX3.5 230v MSB010.CX3.1 110/120v		
Mains L1+N+PE	V/Hz	±10%	230 V 50/60 Hz, 110 V 50/60 Hz
Maximum power consumption	90 W		
Maximum capacity	52.8 ml		
Maximum rotational speed	14500 rpm		
Maximum acceleration	14573 x g		
Maximum kinetic energy	2382 Nm		
Time range (adjustment)	1 - 99 min, or ∞ (1 min. steps)		
Acceleration	3 linear characteristics		
Deceleration	3 linear characteristics		
Programs	9		
PULSE – short duration operation	Yes		
Out-of-balance detection	Yes		
Interference level	PN-EN-55011		
Noise level	56 dB		

Physical data:

Depth	270 mm
Width	220 mm
Height	180 mm
Weight	5 kg

Environmental conditions:

PN-EN-61010-1 (point 1.4.1.)

Ambient temperature	+5°C to +40°C
Relative humidity at ambient temperature	80 %
Installation category	II PN-EN 61010 - 1
Degree of pollution	2 PN-EN 61010 - 1
Protection zone	300 mm

Manufacturer's recommendations

Only the manufacturer-specified centrifuge capillaries, with the correct diameter, length and durability are suitable and should be used for spinning in this centrifuge. The use of equipment made by other manufacturers should be consulted with the manufacturer of the centrifuge. Disinfectants and cleaning agents generally used in medical care should be used in this centrifuge

MSE recommend the use of **Terminex 2** and **Virkon**.

3. INSTALLATION

3.1 Unpacking the centrifuge

Open the package carefully taking care not to damage the machine inside. Take the centrifuge out of the box and ensure that all accessories are checked and identified.

Keep packaging materials at hand for possible transport in the future.

Content of the package (enclosed with every centrifuge).

Cat. No.	Type of accessories	Number Of piece
78500.001	Hexagon key for the rotor	1
76600.020	Emergency release key for opening of the lid	1
87900.001	Rotor Clamp	1
32600.063	Fuses WTA-T 3.15 A 250V	2
35020.151	Power cord 230 V or	1
35020.152	Power cord 110 V	(1)
71100.1219	Operating Instructions	1

3.2 Location

Almost all the energy supplied to the centrifuge is transformed into heat and then emitted to the environment. Therefore, proper ventilation is essential, and ventilation ducts must be kept clear at all times. Moreover, the centrifuge should not be located near any heat source such as radiators and should be kept away from direct sunlight. The bench that the centrifuge sits on should be stable and flat. Around the centrifuge, a 300mm radius safety zone must be established. The normal laboratory operating conditions are from 15°C to 35°C. In the case of changing machine location from extreme cold to warmer conditions, please allow the centrifuge to stand with the lid open for a minimum of 4 hours.

This will reduce any chance of condensation taking place inside the centrifuge.

3.3 Connection to the power supply

Supply voltage has to be consistent with that given on the rating plate. MSE centrifuges come with a 2-metre long IEC lead fitted with plugs resistant to the dynamic loadings.

MSE recommends that an emergency cut -out should be installed far away from the centrifuge, near the exit door from the room or outside the room. Supply voltage should be either 230 V 50/60 Hz, or 110/220 V 50/60 Hz.

3.4 Fuses

The centrifuge has standard protection with the WTA-T 3.15 Amp 250 V fuse, which is situated in the IEC socket at the back of the centrifuge. The fuse fitted in the plug of the IEC lead is 5A.




Figure 1 : Plug-In socket unit

4. SAFETY NOTES

4.1 General information about your centrifuge

The centrifuge has a rigid self-supporting structure. The case is made of Acrylonitrile Butadiene Styrene (ABS) plastic, the front and back is made of aluminum sheet. The rear of the lid is fixed in place by steel hinges and the front of the lid is locked with an electromagnetic lock, which stops the possibility of opening during centrifugation. The centrifuge is an air-cooled machine.

4.2 Operating personnel

The Micro Centaur Plus centrifuge should only be operated by trained laboratory technicians who have read this user manual. The operating manual should be placed near the centrifuge at all times. The centrifuge should not be misused.



4.3 Safety devices

Safety devices described below are incorporated in the design of the Micro Centaur Plus:

4.3.1 Lid lock

The centrifuge will only run when the lid is properly closed. In turn, the lid will only open after the rotor has stopped. In the case of emergency opening of the lid during the operation, the centrifuge is immediately switched-off and the brake is automatically applied until the rotor has completely stopped. With the lid is open, the drive is completely disconnected from the power that makes it impossible to start the centrifuge.

4.3.2 Checking the rotor has stopped

Opening of the centrifuge lid is only possible when the rotor has stopped. This stopped condition is verified by the microprocessor, which sends a signal to the display where the letter "S." appears to show that the rotor has stopped and the lid can be opened. The display will show  when the lid is closed and  when the lid is open.

4.4 Opening the lid manually

It is possible to open the lid manually if there is a power failure or the centrifuge develops a fault. Please use the key provided with your centrifuge to press into the hole on the right-hand side. The lid should now open.

CAUTION: Do not attempt to open the lid whilst the centrifuge is running.

4.5 Good operator practice

Below is a list important points that users should follow in order to get the most out of their centrifuge.

1. Set the centrifuge in horizontal position on a rigid bench.
2. Ensure a safe position for the centrifuge.
3. Ensure free space around the centrifuge (at least 30 cm).
4. Firmly fix the rotor on the motor shaft.
5. Avoid unbalanced loads. If an imbalance has occurred, a “U” will appear next to the acceleration rate.
6. Load the rotor holes with the same tubes as recommended.
7. Do not centrifuge the test tubes of different dimensions.
The possibility to centrifuge test tubes of different dimensions should be avoided at all times. However, if it is absolutely necessary then opposite loads should be the same weight. Care should be taken to load the rotor holes asymmetrically.
8. Ensure that all holes are being used in the rotor.
9. Only use tubes made from a plastic material.
Never use glass test tubes.
11. Fill in all test tubes equally in order to protect the centrifuge against imbalance.
12. Always use tubes that are in good condition and discard any that show signs of wear.
13. Infectious materials should be processed in closed test tubes only.
14. It is prohibited to centrifuge explosive and inflammable materials.
15. It is prohibited to centrifuge substances prone to reaction in result of delivery of high energy during centrifugation.

4.6 Hazards and precautions

1. In order to avoid damaging the centrifuge in any way, the user should read and understand the instruction manual before operating the Micro centaur Plus.
2. Un-trained personnel should never operate the centrifuge.
3. The centrifuge must not be transported with the rotor mounted on the motor shaft.
4. The centrifuge should only be used in conjunction with original rotors and spare parts.
5. If the user discovers a fault with the centrifuge, a qualified service engineer should be contacted through MSE (UK) Ltd.
6. Do not switch the centrifuge on if it is not installed properly or if the rotor is not fitted correctly.
7. The centrifuge must not be operated in places where there is a risk of extreme hot or cold temperatures.
8. Do not use flammable or highly volatile materials in the centrifuge.
9. Users are advised to take extreme care when centrifuging toxic or contaminated substances.
10. The appropriate decontamination procedure should be adhered to when cleaning any parts of the centrifuge that has been in contact with any toxic or hazardous material.
11. Do not open the lid manually when the rotor is still turning.
12. Do not exceed the load limit set by the manufacturer.

13. Rotors are intended for fluids of average homogeneous density equal to 1.2 g/cm³ or smaller when centrifugation is carried out at maximum speed. When fluids of higher density are used, then it may be necessary to reduce the speed.
14. Do not use the rotors that have signs of corrosion or other mechanical defects.
15. Do not use highly corrosive substances, which could cause material impairment and lower the mechanical properties of rotor.
16. Do not use plastic rotors and accessories that have not been recommend by MSE Centrifuges Ltd.
17. Only commercial plastic test tubes should be used.
18. Ensure that the test tubes are in good condition and show no signs of damage or wear.
19. Do not run the centrifuge without the rotor.
20. Do not run the centrifuge without the lid on the rotor.
21. Do not move the centrifuge whilst in use.
22. Maintain a safe distance from the centrifuge whilst in use.
23. Do not place objects on the centrifuge.

4.7 Storage

If your Micro Centaur Plus is stored away for more than 12 months, we recommend that the centrifuge be inspected by an authorised engineer before use.

4.8 Maximum load

In order to avoid overloading of the rotor the operator shall observe the maximum permitted load, which is specified on every rotor. Maximum permissible load is reached when all test- tubes are filled with the fluid with 1.2 g/cm³ density.

If density of the centrifuged liquid is higher than 1.2 g/cm³, then test tubes can only be filled partially. Alternatively, the operator can adjust the speed of the centrifuge by using the following formula.

$$n_{\text{perm}} = n_{\text{max}} * \sqrt{\frac{1,2}{\gamma}}; \quad \gamma = \text{specific gravity } \frac{\left[\frac{\text{g}}{\text{cm}^3} \right]}{\left[\text{cm}^3 \right]}; \quad n_{\text{max}} \text{ [maximum rotational speed - rpm]}$$

5. OPERATING THE CENTRIFUGE

5.1 Getting started

1. Connect the centrifuge to the mains socket.
2. Switch the power supply on by using the on/off switch, situated on the right-hand side of centrifuge.
3. Open the centrifuge lid by pressing the LID key. Prior to fitting the rotor, the user should check that the bowl is free from foreign objects and impurities such as dust and residues of fluid.
4. Mount the rotor on the motor shaft and tighten in a clockwise direction (using the supplied hexagon key).
5. Rotor lids must be fully tightened on the rotor in order to ensure a good seal.
6. Once in use, only use sample containers intended for the rotor.

CAUTION: This Centrifuge compensates for a small weight difference. However, it is recommended that users load the rotor with equal loads to ensure minimal vibrations during operation.

In order to prolong the life of your rotor, lubrication should be used on threaded parts of the rotor and vacuum grease on the “o” rings inside of the rotor lid.

To take the rotor off, loosen the rotor nut by turning in an anti-clockwise direction. Take the rotor off by gently pulling in an upward direction.

5.2 Controls

The microprocessor control system being used in the centrifuge ensures a wide range of parameters and operation can be set:

7. Stepped speed control within the range from 300 to 14000 rpm with increments of 100 rpm or RCF x g.
8. Selection of the centrifuging time within the range from 1 to 99 minutes with 1-minute interval adjustments, or continuous centrifugation.
9. Using the PULSE button to centrifuge for a short period of time.
10. Program up to 9 different operations.
11. 3 different acceleration rates (⏏) to choose from.
12. 3 different de-acceleration rates (⏏) to choose from.

5.2.1 The control panel

The control panel and display is a completely sealed unit that shows key data in a digital format and allows the user to adjust speed, time, acceleration and braking rates. Figure 1.1 shows the layout of the control panel:




Figure 2



The START key is used to start the centrifuge program with the parameters displayed on display. The START key can only be activated when the lid is closed.



The STOP key is used for aborting the operation at any time. After pressing the key an oblique (“▤”) will appear on the display, indicating that the rotor is braking. When the rotor has come to a complete standstill, the letter “S.” will appear on the display signaling that the rotor has stopped. Sequential pressing of the stop key will increase the brake rate to maximum even if the brake rate is set lower.




The LID key is used to open the lid of the centrifuge once it has come to a complete standstill. The LID key can only be used when the display reads “S.” This key cannot be used whilst the rotor is in motion. Once the rotor has stopped the lid can be opened by pressing the lid key, the display will show S▤ and the LED in lid key will be lit up. This denoted that the rotor has stopped and that the lid is open.



The PULSE key is used when the operator wants a short spin without having to program a cycle. The user can hold down this key and the centrifuge will begin to spin. Once the user has let go of the button, the centrifuge will gradually slow down and come to an eventual stop at the brake rate previously chosen. Pressing the STOP key after releasing the PULSE key will increase the brake rate to the maximum even if the set brake rate is lower.

The parameter keys are used to adjust speed, time, RCF, as well as, acceleration and de-acceleration rates. Figure 1.2 shows each of their functions:

Figure 3



Each of the 9 programmes will allow you to set speed, time, acceleration and braking rates.

5.3 Description of sounds

The Micro Centaur Plus will make different 'beep' sounds depending on the information being displayed. Below is a description for each type of beep:

Sound	Description
One short beep	Confirmation that the command has been accepted.
Two short beeps	The command has not been accepted. (e.g. setting the speed above 14 500 rpm)
One long beep	<ol style="list-style-type: none"> 1. Signals the start of braking after pressing the STOP key. 2. Signals the use of the PULSE mode. 3. braking after the use of PULSE mode.
Five short beeps	The rotor has stopped and the lid can be opened
Five short beeps and one long beep	Signals that the centrifuge is ready after switching the power on.

5.4 LCD display

The centrifuge has an LCD display, which shows information regarding the time, speed and RCF, as well as the acceleration and deceleration rates.


When the centrifuge is first switched on, the digital display shows the model name and our web address for three seconds before reverting to the main display.

Figure 4



The main display shows information about the centrifuge cycle. An example is shown below in Figure 5.

Figure 5




^ 1 = acceleration rate 1 (maximum acceleration)

∇ 1 = braking rate 1 (maximum braking rate)

Users can switch between displaying the speed and RCF during the run by pressing the SPEED RCF key twice (see Figure 6)

Figure 6



Note: The operator cannot change the parameters once the centrifuge is in motion.

5.2 Switching the centrifuge on

Having familiarized your self with the instruction manual, you can set the programmers that you require load the rotor, shut the lid and press the START key.

5.3 Switching the centrifuge off

The centrifuge will automatically stop when the programme is completed. It is possible to stop the machine before the programme is completed by pressing the stop key.

When you have completed your requirements, please remember to switch the machine off at the mains switch located on the side of the centrifuge.

6. MAINTENANCE

Please use safety gloves for the following work

6.1 General cleaning of the centrifuge

Before cleaning the centrifuge, users are advised to wear safety gloves to minimize the risk of irritation. For general cleaning, water mixed with a mild detergent is sufficient. Do not use any corrosive substances or abrasive cloths to clean the centrifuge. Do not use alkaline solutions, flammable solvents or agents containing abrasive particles.

6.2 Cleaning of the accessories

It is recommended that rotors are regularly cleaned to ensure optimum performance and prolong their life. This should be done with a non-corrosive cleaning agent (between 6 and 8pH) as these can weaken the rotor. When cleaning the rotor please pay particular attention to the holes of the rotor.

If cracks appear in the rotor, we recommended that the rotor be replaced immediately. Rotors have to withstand high stresses of centrifugal force and therefore any slight damage would be made worse during operation and inadvertently cause damage to the machine.

All accessories used in the centrifuge, such as sample containers, should be clean and free of debris. The sample tubes should also be dry before use in the centrifuge.

6.3 Sterilization and disinfection of the bowl and accessories

When sterilizing the bowl and rotors, users should take into account the material of these items and be sure that the chemicals in the cleaning agents will not react with them. Furthermore, when autoclaving the rotors, users should consider the temperature resistance of the material used in construction. The optimum temperature for autoclaving rotors used in this centrifuge is between 121°C and 124°C. When autoclaving users should also make sure not to exceed 215 kPa (2.15 bar) of pressure. The recommended time for autoclaving is 15 minutes.

The user is responsible for proper disinfections of the centrifuge, particularly if some dangerous material is spilled inside or outside of the centrifuge. During above mentioned works you must always wear safety gloves.

MSE appointed engineers or workshop personnel should not touch the centrifuge unless it has been decontaminated and a certificate is signed by some one in authority to that effect.

6.4 Tube cracking

It is recommended that users operate in a clean environment to reduce the risk of any debris entering your centrifuge system. Debris can crack tubes and rotors and also damage the centrifuge itself.

7. SERVICING YOUR CENTRIFUGE

7.1 Fault finding

Below is a list of faults and their possible solutions:

Fault	Solution
The centrifuge has no display	<ol style="list-style-type: none">1. Check that the power cable is properly plugged in and the mains is turned on2. Check that the fuse hasn't blown in the plug3. Ensure that the power switch is on
The centrifuge does not start	<ol style="list-style-type: none">1. Check that the lid is closed properly2. Ensure that the display is not flashing. The centrifuge could be waiting for you to input some parameters
The centrifuge starts, but it is not accelerating properly	Make sure that the rotor is not overloaded (the letter E will signify this on the display). Turn centrifuge off and on again to clear the error message.
Lid is not opening	The lid lock or lid spring could have failed. Contact MSE for repair.

If any problems persist, please contact your local MSE service agent

7.2 Opening the lid manually

It is possible to open the lid manually if there is a power failure or a fault develops with the centrifuge. Please use the key provided with your centrifuge to press into the hole on the right hand side. The lid should now open.

CAUTION: Do not attempt to open the lid whilst the centrifuge is running

7.3 Safety inspections

From the point of view of operational safety the centrifuge has to be subjected to inspection carried out by authorised service engineer or especially trained experts at least once every year in the state of operational readiness. More frequent inspections could be required if the centrifuge or rotors have been subjected to a corrosion inducing environment. Results of inspections, repairs and tests have to be recorded and kept on file. This instruction manual should be stored near the centrifuge for reference.

7.4 Notes for engineers

The Engineer will pay special attention to the most important parts of the machine and check that they are not damaged in any way. These include:

1. Motor anti vibration ring.
2. The motor shaft.
3. All Centrifuge accessories and especially their structural components. Check they are all free from corrosion.
4. All metal parts. Check they are free from cracks and abrasions.
5. Screws. Check that all screws are tight.
6. Earth connections. Check they are all connected inside the equipment.
7. Rubber sealing rings between the rotor and lid. Check they are in a good condition.

IMPORTANT! Please carry out a safety test for earth bonding, insulation and earth leakage before completing the service.

5.4 Repairing your centrifuge

Your Micro Centaur Plus should provide you with years of trouble-free centrifugation. However, there are occasions when, for various reasons, your centrifuge may need to be repaired.

Warranty repairs will be carried out by authorised service personnel or in an authorised service workshop designated by MSE SAS.

The user will be responsible for the disinfection and decontamination of the unit prior to any service engineer visit or if unit is to be sent back to an MSE Centrifuges Ltd authorised workshop.

The centrifuge will only be accepted if the disinfection / de-contamination certificate is forwarded to the repair shop and a copy of the certificate accompanies the Centrifuge.

Information about authorised service workshops can be obtained from the Manufacturer.

8. HOW TO OBTAIN SERVICE ON YOUR MICRO CENTAUR PLUS

8.1 Warranty

Manufacturer grants the Buyer the guarantee on the conditions specified in the Guarantee Certificate. The Buyer forfeits the right to guarantee repair if using the device inconsistently with the User Manual provisions, when damage results from the User's fault.

Repairs should be carried out in authorized service workshops, granted with the MSE Certificate.

The centrifuge shall be sent to repair after decontaminating disinfections. Information about authorized service workshops can be obtained from the Manufacturer. The Micro Centaur Plus comes with a 2-year warranty from the date of purchase. After the 2 years, an annual service of the centrifuge is recommended. This should be only be carried out by an authorised MSE service engineer. The service life of the centrifuge is approximately 10 years.

8.2 Manufacturer's info

MSE SAS
4B rue du vieux bourg
49340 NUAILLE
FRANCE
tel. (+33) 2 41 64 73 41 (sales department - FRANCE)
e-mail info@mse-sas.fr
website: www.mse-sas.com

If you are located outside the UK, please contact your local distributor. If you do not know who your local distributor is, please contact our sales team via the details above.



MICRO CENTAUR P MICRO CENTRIFUGE

9. Accessories list

MSE Cat No.	Description	Max. RPM	Max. RCF	Max. Capacity	Dimension	Unit
MICRO CENTAUR P						
MSB010.CX3.1	Micro Centaur P laboratory centrifuge (110V/60Hz)	14,500				Each
MSB010.CX3.5	Micro Centaur P laboratory centrifuge (230V/50Hz)	14,500				Each
	<u>ACCESSORIES</u>					
MSE2256	Hematocrite rotor for 24 capillaries	14,500	14,574	24 x 50 mm		Each
MSE6201	Hematocrite reader - round					Each
MSE1252	Angle rotor for PCR strips with Hermetically Sealed Lid (angle 45°)	14,500	13,163	(2 x 8, 2 x 4) x 0.2 ml		Each
MSE1253	Angle rotor 12 x 2/1.5 ml with Hermetically Sealed Lid	14,500	15,279	12 x 2/1.5 ml		Each
MSE4135	Round Adaptor without bottom for 0.5 ml tube				1 x 8 mm	Each
MSE4184	Round Adaptor for round-bottom 0.2 ml tube				1 x 6.2 mm	Each
MSE1255	Angle rotor 24 x 2/1.5 ml with Hermetically Sealed Lid (angle 45°)	14,500	15,279	24 x 2/1.5 ml		Each
MSE4135	Round Adaptor without bottom for 0.5 ml tube				1 x 8 mm	Each
MSE4184	Round Adaptor for round-bottom 0.2 ml tube				1 x 6.2 mm	Each
MSE1254	Angle rotor 18 x 0.5 ml with Hermetically Sealed Lid (angle 45°)	14,000	12,928	18 x 0.5 ml		Each

✂ If required this page may be removed and kept separately.

MICRO CENTAUR PLUS SERVICE REQUEST FORM

Should it become necessary to have your MSE centrifuge repaired. Please take a few moments to fill out this form that will help us to ensure you receive the best and fastest service possible.

Serial number:
(on plate at back of unit)

Date purchased:

Where purchased:

Brief description of fault:
.....
.....

Date fault first occurred:

Date repair centre contacted:

Authorisation number:

Condition of centrifuge:

Has it been disinfected? Yes / No

Disinfectant used:

Contact name:

Address:
.....
.....

Telephone Number:

EU DECLARATION OF CONFORMITY

Medical Scientific Equipment SAS



Declaration Number: **207241-001**

The manufacturer: MSE SAS
Rue du Vieux Bourg
49340 NUAILLE, France

Declare that the product: Laboratory centrifuge

Product name: **Microcentaur + / Microcentaur plus**
Model numbers: **MSB010.CX3.1, MSB010.CX3.5**

Is conformal to the following directives and standards:

Directive 98/79/EC – In Vitro Diagnostic – *Not classified to list A or B and not for self-testing*

Harmonized standards:

EN 15223-1:2016	EN ISO 18113-3:2011
EN 13612:2002	EN 61326-2-6:2006
EN 13612:2002/AC:2002	EN 61010-2-101:2002
EN 13975:2003	EN 62304:2006
EN ISO 14971:2012	EN 62304:2006/AC:2008
EN ISO 1811-3-1:2011	EN 62366:2008

Directive 2014/35/EU – Low voltage (LVD)

Harmonized standards:

EN 61010-1:2010	EN 61010-2-020:2006
-----------------	---------------------

Directive 2014/30/EU – Electromagnetic compatibility (EMC)

Directive 2011/65/EU – RoHS2

This declaration of conformity is issued under our sole responsibility.

Signed for and on behalf of: **MSE SAS**


Place & date of issue: NUAILLE, **30th September 2020**

Name & Position: M. Ricard (**Quality Manager**)



Nomograph

To determine the **RCF (Relative Centrifugal Force)**



DISTRIBUTOR:



MSE SAS

4B rue du vieux Bourg
49340 NUAILLE
FRANCE
Tel: +33 (0) 2 41 64 73 41
E-mail: info@mse-sas.fr
Web: www.mse-sas.com
